

# HLA typing with build 38 and OptiType

Brad Chapman, Miika Ahdesmaki, Justin Johnson  
AstraZeneca Translational Oncology  
Bioinformatics Core, Harvard Chan School

10 February 2016

- 1000 genomes: build 38 + IMGT/HLA-3.18.0
- bwa mem extracts HLA reads
- Map reads only to HLA sequences
- OptiType: Call HLA types

<https://github.com/lh3/bwa/blob/master/README-alt.md#hla-typing>


<https://github.com/FRED-2/OptiType>

<https://github.com/chapmanb/bcbio-nextgen>

# GRCh38 – graph based, many more alternative loci

## Excitement about GRCh38


Alt loci



## Model Centromere Sequences



DPYD  
GGAACGCAG  
GGAACACAG  
R->C



<http://www.slideshare.net/GenomeRef/transitioning-to-grch38>

[ftp://ftp.1000genomes.ebi.ac.uk/vol1/ftp/technical/reference/GRCh38\\_reference\\_genome/](ftp://ftp.1000genomes.ebi.ac.uk/vol1/ftp/technical/reference/GRCh38_reference_genome/)

# Alignment: bwa alternative allele support

Read: ATCAGCATC

```
ALT ctg 1:      TGAAA---CGAATGCAAATGGTCAATCAGCATCGAACTAGTCACAT
                ||||| (high div) ||||| (novel ins) |||||
Chromosome: GCGTACATGATACGAATCgGCATCATGGTC-----CTAGTCACATCGTAATC
                ||||| ||||| (novel ins) |||||
ALT ctg 2:      TGATACGAATCgcCATCATGGTCAATCgcAgCGAACTAGTCACAT
```

4 potential hits: **ATCAGCATC** > **ATCgGCATC** > **ATCgcCATC** > **ATCgcAgC**

2 hit groups: {**ATCAGCATC**, **ATCgcAgC**} and {**ATCgGCATC**, **ATCgcCATC**}

Hits considered in mapQ: **ATCAGCATC** and **ATCgGCATC** (best from each group)

In the output SAM: **ATCgGCATC** as the primary SAM line with mapQ=0

**ATCAGCATC** as a supplementary line with mapQ>0

**ATCgcAgC** as a supplementary line with mapQ>0

**ATCgcCATC** in an XA tag, not as a separate line

<https://github.com/lh3/bwa/blob/master/README-alt.md>

- Map reads to HLA exome 2 and 3 from IMGT
- Matrix of sequence matches to alleles
- Formulate as integer linear program (ILP)
- Use ILP solver, like GNU Linear Programming Kit (GLPK)

<https://github.com/FRED-2/OptiType>

[http:](http://)


[//bioinformatics.oxfordjournals.org/content/30/23/3310](http://bioinformatics.oxfordjournals.org/content/30/23/3310)

- Omixon example data
- Exome (1000 genomes) and deep targeted data
- HLA type I calls (A, B, C)
- Good validation results
  - 24/24 (100%) on targeted
  - 22/24 (92%) on exome

<http://www.omixon.com/hla-typing-example-data/>

<https://gist.github.com/chapmanb/8f994618a7fc5e88f893>

# HLA P-group resolution



© SGE Marsh 04/10

<https://www.ebi.ac.uk/ipd/imgt/hla/>

[http://hla.alleles.org/alleles/p\\_groups.html](http://hla.alleles.org/alleles/p_groups.html)


- Build 37 and 38
- Validation sets: Genome in a Bottle, Illumina Platinum Genomes
- Lift-over methods: CrossMap/LiftOver, NCBI Remap
- 38 builds: with/without alternative alleles
- Variant callers: FreeBayes, GATK  
HaplotypeCaller

<http://bcb.io/2015/09/17/hg38-validation/>




# GRCh38 – advantage for variant calling

## Reference assembly influence




# Avoiding collapsed repeats

Distribution of variants on chromosome 1



Distribution of variants on chromosome 1



<http://www.slideshare.net/kmsteinberg/>

the-importance-of-high-quality-reference-genome-assemblies-to-personal-and-medical-genomics



Genome in a Bottle  
Consortium



**Global Alliance**  
for Genomics & Health

ICGC-TCGA DREAM Mutation Calling challenge

<http://www.genomeinabottle.org/>

<http://ga4gh.org/#/benchmarking-team>

<https://www.synapse.org/#!Synapse:syn312572>

# hg19/hg38 comparison: NA12878 Platinum Genomes

## SNPs: freebayes



## SNPs: gatk-haplotype



0% 0.2% 0.4% 0.6% 0.8% 1% 1.2% 1.4%

0% 0.2% 0.4% 0.6% 0.8% 1% 1.2% 1.4%

## Indels: freebayes



## Indels: gatk-haplotype



0% 2% 4% 6% 8% 10% 12%


0% 2% 4% 6% 8% 10% 12%

False negative rate


False discovery rate

# GRCh37/hg38 comparison: NA12878 Genome in a Bottle


## SNPs: freebayes



## SNPs: gatk-haplotype



## Indels: freebayes



## Indels: gatk-haplotype



False negative rate

False discovery rate

- SNPs: build 38 more sensitive
- SNPs: build 38 reduces false positives
- Indels: build 38 detected more
- Indels: work on sensitivity and precision

Need conversion approaches for resources not yet available on build 38

- CrossMap:

<http://crossmap.sourceforge.net/>

- NCBI remap:

<http://www.ncbi.nlm.nih.gov/genome/tools/remap>

- Both performed well

- NCBI remap has additional sensitivity, but requires tuning